Extra Wide Elevated Leg Rest

A good day starts with a great night's sleep. Back pain, knee pain, varicose veins, restless leg syndrome, and other medical conditions can often rob you of the sleep your body needs. The Kolbs Extra Wide Elevated Leg Rest offers optimal positioning for a comfortable and relaxing slumber.

The Kolbs Extra Wide Elevated Leg Rest elevates your legs comfortably above your heart for better blood circulation, and less pressure on your back. Sleep healthy, sleep well.

Medical Grade Foam Core

A high density core layer to ensure just the right firmness and proper support.

Luxurious Jacquard Pillow Case

A removable, jacquard pillow case protects the Elevated Leg Pillow from stains while promoting a comfortable night's sleep.

Foam Wedge is Ideal for
Foot or Leg Elevation

Care Instructions: Cover can be removed and machine washed on delicate cycle. Maximum temperature 60°C (140°F). Do not bleach. Do not iron. Do not dry clean. Tumble dry low or air dry. Foam core can be spot cleaned with a moist, soapy rag. Let air dry. Foam core may never be placed in the washing machine.

Important Notice: The Extra Wide Elevated Leg Rest is shipped in compressed packaging. Remove Cushion from plastic packaging to begin decompression. The process will begin immediately. Allow several hours for cushion to decompress before use.

Please note: Many factors contribute to the speed of your cushion's decompression. In the vast majority of cases your cushion will decompress in the time frame indicated. In the rare case that it does not, please allow the cushion to decompress overnight.

Please also note: Although you may use the cushion at this point, it may take up to 24-48 hours for the leg pillow to reach its peak firmness.

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